

Works: - E-41, ROAD No.2, SARDAR INDUSTRIAL ESTATE, AJWA ROAD, BARODA 19 GUJARAT. INDIA
Office: - 2 TILAK NAGAR SOCIETY, OPP. NAVJIVAN, AJWA ROAD, BARODA 19 GUJARAT. INDIA
Telefax: - 0265 2511484, **(m)** 9979 884940, **(e)** info@bororlab.in, **(web)** www.bororlab.in

ROTARY EVAPORATOR (LABORATORY MODEL)



Motor & Gear Box
Rotation Speed

: Durable Gear Box with brushless motor for minimum noise
: Digital RPM Display & VFD based Speed Control, 0-98 R.P.M

Heating Bath

: SS304 Jacketed Bath, with Digital Temperature Indicator & Controller.

Temperature Range
of Heating Bath

: 0 to 130 Deg. C
: Automatic Motorized VFD based Up & Down motion of Heating Bath.

Condenser

: All SS304 (pharma grade) body.

Power Supply

: Glass Condenser with 0.010m² H.T.A.

: 230 V/ 50 Hz

SAILABLE FEATURES

- Comes with LCD screen.
- All glass fittings in Laboratory Joints (cone & socket)
- Auto controlled digital display of rotational speed (RPM) & bath temperature.
- Digital display of process time.
- Automatic bath lifting.
- Withstands full Vacuum.
- Ideally suited for heat sensitive material.